


NEW



TUBULEX®

- Hastens the healing process after extraction.
- Creates an antiseptic zone of approx. 30 - 40 mm.
- Minimizes pain and inflammatory symptoms.
- Prevents the emergence of "dry socket syndrome".
- Used to disinfect the mucous membrane in the mouth.
- Disposable sterile dressing.



1. Open the aluminium foil plate by pressing with the thumb.
2. Take out a dressing using tweezers.
3. Fold the dressing into a "V" shape before applying.
4. Apply 1-2 dressings to cover the alveolar.